# Decimation of EM Volume Mesh Objects with Pyvista

Examples from MICrONS and H01 Volumes

#### Volume Examples

pinky

MICrONS Layer 2/3 mouse visual cortex

minnie

■ MICrONS 1mm³ mouse visual cortex

H01

Harvard/Google human temporal cortex

#### **Decimation Process**

#### Download

 Download mesh files to local disk (or Google Drive for Google Colab)

#### Decimate

 Process and decimate mesh files (creates a new ply file 95% smaller)

#### Visualize

 Visualize decimated mesh objects (all examples shown here use vtk)


#### Source: downloading meshes

- Mesh download code from Forrest Collman's <u>MeshExample.ipynb</u>
- Thank you Bethanny Danskin (Allen Institute) for providing a workaround for meshparty legacy file handling error
- See <u>Allen Institute Microns Binder github repository</u> for more details and additional notebooks for working with the MICrONS volume datasets

#### Source: pyvista mesh decimation

- Mesh decimation code from Tyler Sloan's <u>Mesh</u> <u>Decimation Pipeline.ipynb</u>
- See <u>Quorumetrix github repository</u> for more details and additional scripts for working with volume meshes in Blender (however, all visualization examples shown here use vtk/OpenGL)

#### MICrONS Layer 2/3 (pinky) volume



Explore EM data from layer 2/3 of the mouse visual cortex

#### Mesh Types in Layer 2/3 Volume

Cell

Cell soma meshes

Mito

Mitochondria meshes

Nuc

Nuclei meshes

#### Layer 2/3: download meshes



Visit the <u>download\_meshes</u> folder on github



Use the respective notebooks to download and view the cell, mitochondria, and nuclei meshes



The cell body and nuclei notebooks also visualize the downloaded meshes with matplotlib and plotly plots

```
# setup the mesh meta to handle downloads and
caching
mesh dir = 'data/neuron meshes v185/'
seg_source =
"precomputed://gs://microns public datasets/pinky100
v185/seg"
cv obj = cloudvolume.CloudVolume(seg source,
use https=True)
seg id = [648518346349525821]
```

Workaround code for meshparty legacy issue with Layer 2/3 cell body mesh files


```
# Iterating through the seg_id list
for seg in seg_id:
    print(f"Processing seg_id: {seg}")

# Retrieve the mesh object with tqdm progress bar
    print("Retrieving mesh...")
    cv_mesh = cv_obj.mesh.get(seg,
remove_duplicate_vertices=True)
```

```
# Extracting and reshaping faces if necessary
print("Processing mesh faces...")
faces = np.array(cv_mesh.faces)
if len(faces.shape) == 1:
    faces = faces.reshape(-1, 3)
# Creating the Mesh object
print("Creating Mesh object...")
mesh = Mesh(vertices=cv mesh.vertices, faces=faces)
```

```
# Constructing the file path
mesh file = os.path.join(mesh dir, str(seg) + '.h5')
# Writing the mesh to file
print(f"Writing mesh to file: {mesh file}")
mesh.write to file(mesh file)
print(f"Finished processing seg id: {seg}\n")
```

matplotlib visualization of a pyramidal cell mesh



#### Mesh download for mitochondria

# setup the mesh meta to handle downloads and caching

```
mito_mesh_dir = 'data/meshes/'
mito_source =
"precomputed://https://td.princeton.edu/sseung-archive/pinky100-mito/seg_191220"
mito_mm =
trimesh_io.MeshMeta(cv_path=mito_source,
disk_cache_path=mito_mesh_dir)
```

#### Mesh download for mitochondria

```
cellid_list = [648518346349529031]

mitodf =
mito[mito['cellid'].isin(cellid_list)]
mitodf_mitoid_list = mitodf.mito_id.to_list()
```

#### Mesh download for mitochondria

```
for i in range(len(mitodf mitoid list)):
    mito id = mitodf mitoid list[i]
    mito seg id = mito id
    mito downloadmesh = mito_mm.mesh(seg_id =
mito_seg_id, remove duplicate vertices=True)
    print(f'completed ' + str(i) + ' of ' +
str(len(mitodf mitoid list)))
```

#### Mesh download for nuclei


```
nuc_mesh_dir = 'data/nuclei meshes/'
nuc source =
"precomputed://https://td.princeton.edu/sseung
-archive/pinky100-nuclei/seg"
nuc mm =
trimesh io.MeshMeta(cv path=nuc source,
disk cache path=nuc mesh dir)
```

#### Mesh download for nuclei

```
nuc_id = [5645]

for i in tqdm(range(len(nuc_id)),
  desc="Downloading Meshes"):
  downloadmesh = nuc_mm.mesh(seg_id=nuc_id[i],
  remove_duplicate_vertices=True)
```

plotly visualization of a nucleus mesh



#### Layer 2/3: decimate meshes



Visit the <u>decimated\_meshes</u> folder on github




Use the respective notebooks to decimate the cell, mitochondria, and nuclei meshes




Options for single and batch processing meshes using lists

#### Mesh decimation: nucleus example

Original Mesh Decimated Mesh






#### Layer 2/3 volume: visualize meshes




Visualize decimated meshes using a variety of methods with vtk with customized notebooks in the <u>decimated\_meshes</u> directory




Additional options include adding pre- and post-synaptic sites, filtered for pre-synaptic sites for a neuron of interest, highlighted cells, highlighted mitochondria by size, and more



All mitochondria in 44 astrocytes in the Layer 2/3 volume with vasculature shown in red



#### MICrONS 1 mm<sup>3</sup> (minnie) volume



#### minnie volume: download meshes



Visit the Jupyter Notebook folder on github




Use <u>minnie\_plotly\_and\_vtk\_visualizer.ipynb</u> to download and view cell meshes (list-based approach)

```
# setup the mesh meta to handle downloads and
caching
mesh dir = 'data/cell meshes'
seg source =
"precomputed://gs://iarpa microns/minnie/minnie65/
seg"
mm = trimesh io.MeshMeta(cv path=seg source,
                         disk_cache_path=mesh dir,
                         cache size=20)
```

```
cell id list = [864691134964446239,
864691135012524790,
                    864691135256138671,
864691135269913253,
                    864691135337851366,
864691135341421745,
                    864691135348239831,
864691135386363265,
                    864691135497750291,
864691135657783170,
                    864691135809446092,
864691135837182867,
                    864691135880405261,
864691136023933241,
                    864691136194008022,
864691136330394007
```

```
for i in tqdm(range(len(cell_id_list)),
desc="Downloading Meshes"):
    downloadmesh = mm.mesh(seg_id=cell_id_list[i],
    remove_duplicate_vertices=True)
```

### Plotly visualization of mesh objects



#### minnie volume: decimate meshes



Decimate meshes using a list of cellids of interest with <a href="mailto:pyvista\_decimate\_cell\_bodies\_minnie.ipynb">pyvista\_decimate\_cell\_bodies\_minnie.ipynb</a>

#### minnie volume: visualize meshes




Visualize <u>cell meshes</u> by cell type (e.g., neurons, astrocytes, vasculature) with vtk

vtk OpenGL rendered neurons and astrocytes

#### Harvard/Google H01 volume

H01 Release Explore Gallery Manuscript Layers SegCLR Embeddings CREST Proofreading Released Data Code



#### H01 volume: download meshes



Visit the Jupyter Notebook folder on github



Start with <a href="https://www.neshes">h01\_plotly\_visualizer.ipynb</a> to download and view cell meshes




See <a href="https://holder.com/holler.gov/holler.com/h

```
# setup the mesh meta to handle downloads and
caching
mesh dir = 'data/c3 cell meshes/' # or change to
your desired folder
seg source = "precomputed://gs://h01-
release/data/20210601/c3"
mm = trimesh io.MeshMeta(cv path=seg source,
                         disk cache path=mesh dir,
                         cache size=20)
```

```
cell_id_list = [28541451958,27843438291,
4970736617]

for i in tqdm(range(len(cell_id_list)),
desc="Downloading Meshes"):
    downloadmesh = mm.mesh(seg_id=cell_id_list[i],
remove_duplicate_vertices=True)
```

## Plotly visualization of mesh objects



#### H01 volume: decimate meshes



Decimate <u>astrocyte meshes</u> of interest as a list (can be used for other cell types too)



Complex meshes, such as astrocytes and vasculature, may take hours (~12 hours per astrocytes) or run out of memory on a typical computer



Use <u>Google Colab</u> with High-RAM (subscription to Google Colab Pro required) to allocate enough memory and reduce processing time (<10 min)

#### H01 volume: visualize meshes





Visualize <u>decimated astrocyte meshes</u> using a gradient coloring method with vtk



Visualize <u>decimated meshes by four cell categories</u> (neurons, astros, microglia, and unidentified segments) with vtk

Gradient coloring of astrocytes





#### Acknowledgements: MICrONS



Read the original research papers on the <u>Citation</u> page at Allen Institute



Read the Terms and Conditions page



Use under Creative Commons by Attribution 4.0 International

#### Acknowledgements: H01



Link to the original research paper on the <u>Manuscript</u> page of H01 Release website



Read the license on the Released Data page



Use under Creative Commons by Attribution 4.0 International

#### Code Availability

Visit my github repositories for each volume:

- layer23-volume
- minnie-volume
- h01-volume